FREEPRINT® MODEL PRO

DETAX
 Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
 ☎ +49 7243/510-0 · 📠 +49 7243/510-100
 www.detax.com · post@detax.com

Made in
 Germany
 04/2024
 5

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FREEPRINT® MODEL PRO

ENGLISH

Light-curing resin based on (meth)acrylate, for the generative fabrication of dental models for DLP printers with UV-LED 385 nm

Area of application

3D Print working and presentation models for the entire spectrum of models in dentistry

SUITABLE FOR THE FOLLOWING DLP PRINTERS/CLEANING/POST-EXPOSURE

see "Annex 1" (enclosed separately)

Processing

- ▶ The properties of the final product depend, among other things, on post-processing. Therefore it must be ensured that the light unit is in an orderly condition and that the models are completely cured (refer to the process description).
- ▶ Shake the material intensively before use and homogenise with a bottle roller.
- ▶ Processing temperature 23 °C ± 2 °C.

Safety Information

- ▶ Only to be used by trained specialists for the specified area of application.
- ▶ Product contains (meth)acrylates and phosphine oxides. The ingredients of **FREEPRINT® MODEL PRO** may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further.
- ▶ Avoid direct contact with the liquid material and the components before post-curing.
- ▶ Irritating to airways, eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ In case of contact with skin, wash immediately with plenty of soap and water.
- ▶ The processing instructions and precautionary measures must be observed.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

Notes

- ▶ DETAX does not accept liability for any damage caused by misuse.
- ▶ Always keep container tightly closed, carefully close immediately after each use.
- ▶ Product may cause allergic reactions.
- ▶ Read and understand the safety data sheet!

Storage

- ▶ **FREEPRINT® MODEL PRO** is to be stored dry (at 15 °C - 28 °C) and protected from light. Minimal influence of light can already induce polymerisation.
- ▶ To protect against contamination, cover the material in the tray with the lid or a glass plate.

Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Manufacturing process

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

Construction process

Generation of a Print Job complying with machine and material parameters

Post-processing


After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

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