
FREEPRINT[®]

3D PRINTING MATERIAL

DENTAL MED RESIN

» **CROWN** «

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Composite for dental 3D printing

INDICATION

Permanent single crowns (anterior and posterior areas), permanent inlays and onlays, veneers, long-term provisional crowns, bridges and permanent denture teeth.

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see "Annex 1" (enclosed separately)

PROCESSING

- ▶ The properties of the final product depend, among other things, on the finishing process. Correct post-exposure is important for biocompatibility. Therefore, it must be ensured that the exposure unit is in proper condition and the molded parts are completely cured (refer to the process description).
- ▶ Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- ▶ Do not use thermal methods for disinfection or sterilization. This might deform the workpiece.
- ▶ **FREEPRINT® CROWN** is used for the fabrication of permanent single crowns in the anterior and posterior region, denture teeth, inlays and onlays. It is also used for long-term provisional crowns and bridges having up to 3 pontics with an intermediate pontic and veneers. The minimum occlusal thickness is 1.5 mm, circularly 1.5 mm and cervically 1.0 mm. The connection cross-section for intermediate pontics of a bridge must be at least 16 mm². Failure to observe this information may lead to an undesirable result.
- ▶ A layer thickness of 50 µm is suitable for the printing process.
- ▶ Mechanically polish the surface. Pre-polishing is done with rotating brushes and pre-polishing paste, a high-gloss polish with buffing and high-gloss polishing agents for resin composite.
- ▶ The restorations can be used with the commonly used permanent luting cements, e.g. Variolink Esthetic DC by Ivoclar or RelyX Unicem2 by 3M Espe. Observe the manufacturer's instructions.
- ▶ **Attachment of permanent denture teeth**
Once the denture base (see additional IFU **FREEPRINT® DENTURE**) and the denture teeth have been removed from the respective build platforms and the corresponding cleaning processes completed, leave the printed components in their green body state (not post-exposed). Place the printed denture teeth in the corresponding tooth receptacles of the denture base, check the correct seating of the teeth and remove them again. Apply liquid **FREEPRINT® DENTURE** to the receptacles with a small brush or a syringe, and connect the repositioned denture teeth by exposure in a xenon flash unit (Otoflash G171) with 400 flashes under inert gas conditions (nitrogen). Finally post-cure the complete denture in a xenon flash unit (Otoflash G171) with 2 × 1800 flashes under inert gas conditions (nitrogen).
- ▶ Processing temperature 23 °C ± 2 °C.

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The Summary Safety and Clinical Performance Report (SSCP) for BASIS UDI-DI ++EDET0996K (FREEPRINT® CROWN) is stored in the European Database for Medical Devices (EUDAMED) under the internet link <https://ec.europa.eu/tools/eudamed>.

Until full functionality of the EUDAMED is announced, the SSCP can be requested from regulatory.affairs@detax.com.

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Finishing process

After raising the platform, a dripping time of approx. 10 minutes is recommended. If possible, the finishing should take place immediately after the building process.

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Composite for dental 3D printing

INDICATION

Single crowns (anterior and posterior areas), inlays and onlays, veneers for a wearing time up to 5 years, long-term provisional crowns, bridges for a wearing time up to 12 month and permanent denture teeth.

PATIENT TARGET GROUP

Persons being treated in the context of a dental procedure.

INTENDED USERS

Dentist, dental technician

SUITABLE FOR THE FOLLOWING DLP PRINTERS/CLEANING/POST-EXPOSURE

see "Annex 1" (enclosed separately)

PROCESSING

- ▶ The properties of the final product depend, among other things, on the finishing process. Correct post-exposure is important for biocompatibility. Therefore, it must be ensured that the exposure unit is in proper condition and the molded parts are completely cured (refer to the process description).
- ▶ Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- ▶ Do not use thermal methods for disinfection or sterilization. This might deform the workpiece.
- ▶ **FREEPRINT® CROWN** is used for the fabrication of permanent single crowns in the anterior and posterior region, denture teeth, inlays and onlays. It is also used for long-term provisional crowns and bridges having up to 3 units with an intermediate pontic and veneers. The minimum occlusal thickness is 1.5 mm, circularly 1.5 mm and cervically 1.0 mm. The connection cross-section for intermediate pontics of a bridge must be at least 16 mm². Failure to observe this information may lead to an undesirable result.
- ▶ A layer thickness of 50 µm is suitable for the printing process.
- ▶ Mechanically polish the surface. Pre-polishing is done with rotating brushes and pre-polishing paste, a high-gloss polish with buffing and high-gloss polishing agents for resin composite.
- ▶ The restorations can be used with the commonly used permanent luting cements, e.g. Variolink Esthetic DC by Ivoclar or RelyX Unicem2 by 3M Espe. Observe the manufacturer's instructions.
- ▶ **Attachment of permanent denture teeth**
Once the denture base (see additional IFU **FREEPRINT®DENTURE**) and the denture teeth have been removed from the respective build platforms and the corresponding cleaning processes completed, leave the printed components in their green body state (not post-exposed). Place the printed denture teeth in the corresponding tooth receptacles of the denture base, check the correct seating of the teeth and remove them again. Apply liquid **FREEPRINT® DENTURE** to the receptacles with a small brush or a syringe, and connect the repositioned denture teeth by exposure in a xenon flash unit (Otoflash G171) with 400 flashes under inert gas conditions (nitrogen). Finally post-cure the complete denture in a xenon flash unit (Otoflash G171) with 2 × 1800 flashes under inert gas conditions (nitrogen).
- ▶ Processing temperature 23 °C ± 2 °C.

SAFETY INFORMATION

- ▶ Only for the specified intended use by trained specialists.
- ▶ Avoid direct contact with the liquid material and the components before post-curing. Irritating to airways, eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (e.g. protective gloves, goggles) when handling the uncured material.

- ▶ For finishing the cured material, suitable personal protective equipment must be used (protective gloves, goggles, face mask).
- ▶ In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ In case of contact with skin, wash immediately with plenty of soap and water.
- ▶ Biocompatibility is ensured only after complete polymerization.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

NOTES

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly closed, carefully close immediately after each use.
- ▶ Read and understand the safety data sheet!

For users and/or patients:

Any serious incident that has occurred in connection with this product must be reported immediately to incident@detax.com and to the competent public authority of the Member State in which the user and/or patient resides.

Brief Report on Safety and Clinical Performance (SSCP):

The Summary Safety and Clinical Performance Report (SSCP) for BASIS UDI-DI ++EDET0996K (freeprint® crown)

is stored in the European Database for Medical Devices (EUDAMED) under the internet link

<https://ec.europa.eu/tools/eudamed>. Until full functionality of the EUDAMED is announced, the SSCP can be requested from regulatory.affairs@detax.com.

STORAGE

- ▶ Store FREEPRINT® CROWN in a dry place (at 15 °C – 28 °C) and away from light. Even slight exposure to light can trigger polymerization.
- ▶ To protect against contamination, cover the material in the tray with the lid or a glass plate.

CONTRAINDICATION

Contains (meth)acrylates and phosphine oxides.

Ingredients of FREEPRINT® CROWN may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further. Apply FREEPRINT® CROWN intraorally only when fully polymerized.

SIDE EFFECTS

Product may cause allergic reactions.

DISPOSAL

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Storage:

Processing:

at 23 °C ± 2 °C

Manufacturing process

Data processing and creation of the support structure according to the specifications of the CAD software manufacturer

Building process

Creation of a print job in compliance with the machine and material parameters, see table

Finishing process

After raising the platform, a dripping time of approx. 10 minutes is recommended. If possible, the finishing should take place immediately after the building process. The 3D-printing resin can remain in the vat for a max. of 14 days.

Pre-cleaning

Remove components from the platform and remove excess resin with a gentle stream of (pressurized) air.

Cleaning

see "Annex 1, CLEANING EQUIPMENT"

Then clean cavities and surfaces additionally with compressed air and turn the component over.

Remove support structures from the component.

Main cleaning

see "Annex 1, CLEANING EQUIPMENT"

Dry cavities and surfaces with compressed air. The surface must no longer be shiny, otherwise resin residues will be present.

Drying

Store the component protected from light at room temperature for 30 minutes before post-exposure.

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*The set of parameters includes all relevant material and printer specific information

CLEANING EQUIPMENT

Cleaning unit Manufacturer, Model	Cleaning process
Ivoclar PrograPrint Clean	Clean the parts with isopropyl alcohol (purity $\geq 98\%$) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 2 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Rapidshape RS wash	Use the following settings: DETAX Freeprint-crown Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Ultrasonic bath Bandelin Sonorex	Clean the parts with isopropyl alcohol (purity $\geq 98\%$) for 1 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 1 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.

LIGHT CURING EQUIPMENT

Light Curing unit Manufacturer, Model	Curing process
Ivoclar PrograPrint Cure	Post curing A: Wavelength = 405 nm; Intensity = 45 %; Duration = 5 min Post curing B: Wavelength = 460 nm; Intensity = 45 %; Duration = 5 min
NK Optik Otofash G171	2 x 2000 flashes under inert gas, turnaround components after 2000 flashes
Rapidshape RS Cure	Use the following settings: DETAX freeprint-crown

DETAX

DETAX GmbH

Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany

+49 7243/510-0 · +49 7243/510-100

www.detax.com · post@detax.com

Made in
Germany

06/2024

7