FREEPRINT[®]

3D PRINTING MATERIAL

SURGICAL GUIDES

DETAX

DE	Gebrauchsanweisung	3	NO	Bruksanvisning	47
EN	Instructions for use	7	PL	Instrukcja stosowania	51
FR	Mode d'emploi	11	PT	Manual de instruções.....	55
ES	Modo de empleo	15	RO	Instrucțiuni de utilizare	59
IT	Istruzioni per l'uso	19	SV	Bruksanvisning	63
DA	Brugsanvisning.....	23	SK	Návod na použitie.....	67
EL	Οδηγίες χρήσης.....	27	SL	Navodila za uporabo	71
HR	Upute za uporabu	31	CS	Návod k použití	75
LV	Lietošanas instrukcija.....	35	HU	Használati útmutató.....	79
LT	Naudojimo instrukcija	39	RU	Инструкция по применению	83
NL	Gebruiksaanwijzing	43	TR	Kullanım kılavuzu	87

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Resin for dental 3D printing

INDICATION

Orthodontic bases, surgical guides and X-ray templates for implant dentistry, occlusal orthos, fixation and transfer keys

PATIENT TARGET GROUP

Persons being treated in the context of a dental procedure.

INTENDED USERS

Dentist, dental technician

SUITABLE FOR THE FOLLOWING DLP PRINTERS/CLEANING/POST-EXPOSURE

see "Annex 1" (enclosed separately)

PROCESSING

- ▶ The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the moulds are completely cured (observe process description on page 10).
- ▶ Bottle: Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- ▶ Eco Bag: Before use, the material should be shaken intensively and homogenized by tumbling or using a roller mixer with a suitable attachment.
- ▶ Minimum material thickness for the design is 1.5 mm
- ▶ Polish surface mechanically.
- ▶ Processing temperature $23\text{ °C} \pm 2\text{ °C}$.

FREEPRINT® ORTHO 385 - STERILIZATION OF DRILLING TEMPLATES

- ▶ Locally applicable regulations and hygiene guidelines are to be observed.
- ▶ Important: after main cleaning and prior to post-curing, the components to be sterilized are to be dried, e.g. to be heated in an oven to approx. 40 °C for 30 minutes to remove remaining solvent from the cleaning process
- ▶ FREEPRINT® ORTHO 385 is validated for the „W&H Lisa 522, Program UNIVERSAL 121“ sterilization process (technical data on the right). Steam sterilization may only be performed with equipment that complies with the standards EN 13060 and EN 285.
- ▶ Responsibility for sterility lies with the user. Safety Information

FREEPRINT® ORTHO (405 NM)

- ▶ Do not use heat-based methods for disinfection or sterilization. This could possibly deform the workpiece.

TECHNICAL DATA FOR THE STERILISATION PROCESS

- ▶ A fractionated pre-vacuum with a plateau phase at 121 °C/ 15 minutes is to be selected, whereby this refers to the requirements for a B-cycle small steriliser in terms of its venting performance.

SAFETY INFORMATION

- ▶ Only for the specified intended use by trained specialists.
- ▶ Avoid direct contact with the liquid material and the components before post-curing, especially in pregnant / breastfeeding women. Irritating to eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- ▶ After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- ▶ After contact with skin wash immediately with water and soap.
- ▶ Biocompatibility is only guaranteed with complete polymerisation.
- ▶ MMA-free.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

CLEANING

By dental laboratory / practice

In dental laboratories and practices, the cleaning is done in ultrasonic baths.

By the patient

For daily cleaning, use lukewarm water, a soft toothbrush, liquid soap or detergent. Then rinse well afterwards.

Cleaning tablets can be used following the manufacturer's instructions.

Note

- ▶ Do not use toothpaste; the abrasive particles are coarse and scratch the surface.
- ▶ Do not use any acids or solvents.

NOTE

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly sealed, immediately close the container carefully after each use.
- ▶ Please follow the instructions on the safety data sheet!

For users and/or patients:

Any serious incidents occurring in relation to this product should be reported immediately to incident@detax.com and to the competent authority of the Member State in which the user and/or patient is established.

STORAGE

- ▶ Store **FREEPRINT® ORTHO** in a dry place (at 15 °C – 28 °C) and away from light. Even slight exposure to light can trigger polymerization.
- ▶ To protect against contamination, cover the material in the tray with the lid or a glass plate.

CONTRAINDICATION

Contains (meth)acrylics and phosphine oxide.

Some ingredients of **FREEPRINT® ORTHO** may cause allergic reactions in predisposed persons. In such cases refrain from using the product. **FREEPRINT® ORTHO** only insert intraorally in completely polymerised state.

ADVERSE EFFECTS

Product may cause allergic reactions.

DISPOSAL

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Storage:

Application:

At 23 °C ± 2 °C

Manufacturing

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

Construction process

Generation of a Print Job complying with machine and material parameters

Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended. If possible, post-processing should commence immediately following the construction process.

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*The set of parameters includes all relevant material- and printer specific information

CLEANING EQUIPMENT

Cleaning unit Manufacturer, Model	Cleaning process
Rapidshape RS Wash	Use the following settings: DETAX Freeprint-ortho uv Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Ultrasonic bath Bandelin Sonorex	Clean the parts with isopropyl alcohol (purity \geq 98 %) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity \geq 98 %) for 3 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.

LIGHT CURING EQUIPMENT

Light Curing unit Manufacturer, Model	Curing process
NK Optik Otoflash G171	2 x 2000 flashes under inert gas, turnaround components after 2000 flashes
NK Optik Otoflash 250/500	4000 flashes under inert gas @15 Hz
Rapidshape RS Cure	Use the following settings: DETAX Freeprint-ortho uv

DETAX

DETAX GmbH
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
☎ +49 7243/510-0 · 📠 +49 7243/510-100
www.detax.com · post@detax.com

Made in
Germany
09/2023
3