FREEPRINT®

3D PRINTING MATERIAL

INDIVIDUAL
FUNCTIONAL TRY-INS

DETAX

DE	Gebrauchsanweisung	3	PT	Manual de instruções.....	42
EN	Instructions for use	6	RO	Instrucțiuni de utilizare	45
FR	Mode d'emploi	9	SV	Bruksanvisning	48
ES	Modo de empleo.....	12	SK	Návod na použitie.....	51
IT	Istruzioni per l'uso.....	15	SL	Navodila za uporabo	54
DA	Bruksanvisning.....	18	CS	Návod k použití	57
EL	Οδηγίες χρήσης.....	21	HU	Használati útmutató.....	60
HR	Upute za uporabu	24	RU	Инструкция по применению	63
LV	Lietošanas instrukcija.....	27	TR	Kullanım kılavuzu	66
LT	Naudojimo instrukcija	30	FI	Käyttöohje.....	69
NL	Gebruiksaanwijzing	33	ET	Kasutusjuhend	72
NO	Bruksanvisning	36	BG	Инструкция за употреба.....	75
PL	Instrukcja stosowania	39			

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Resin for dental 3D printing

INDICATION

Individual functional try-ins of complete dentures

PATIENT TARGET GROUP

Persons being treated in the context of a dental procedure.

INTENDED USERS

Dentist, dental technician

SUITABLE FOR THE FOLLOWING DLP PRINTERS

"Annex 1" (enclosed separately)

PROCESSING

- ▶ The properties of the final product depend, among other things, on the finishing process. Correct post-exposure is important for biocompatibility. Therefore, it must be ensured that the exposure unit is in proper condition and the molded parts are completely cured (refer to the process description).
- ▶ Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- ▶ A layer thickness of 100 µm is suitable for the printing process.
- ▶ Polish the surface mechanically if necessary. Pre-polishing is done with rotating brushes and pre-polishing paste, high-gloss polishing with buffing and high-gloss polishing agents for resin.
- ▶ Processing temperature 23 °C ± 2 °C.

SAFETY INFORMATION

- ▶ Only for the specified intended use by trained specialists.
- ▶ Avoid direct contact with the liquid material and the components before post-curing, in particular in pregnant / breastfeeding women. Irritating to airways, eyes and skin (sensitization possible).
- ▶ Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- ▶ Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- ▶ In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- ▶ In case of contact with skin, wash immediately with plenty of soap and water.
- ▶ Biocompatibility is guaranteed only after complete polymerization.
- ▶ Refer to the relevant safety data sheet for hazard and safety information.

NOTES

- ▶ DETAX shall not be held liable for any damage caused by misuse.
- ▶ Always keep container tightly closed, carefully close immediately after each use.
- ▶ Read and understand the safety data sheet!

For users and/or patients:

Any serious incident that has occurred in connection with this product must be reported immediately to incident@detax.com and to the competent public authority of the Member State in which the user and/or patient resides.

STORAGE

- ▶ Store **FREEPRINT® tryin** in a dry place (at 15 °C – 28 °C) and away from light. Even slight exposure to light can trigger polymerization.
- ▶ To protect against contamination, cover the material in the tray with the lid or a glass plate.

CONTRAINDICATION

Contains (meth)acrylates and phosphine oxides.

Ingredients of **FREEPRINT® tryin** may cause allergic reactions in accordingly disposed persons. In such a case, do not use the product any further. Apply **FREEPRINT® tryin** intraorally only when fully polymerized.

SIDE EFFECTS

Product may cause allergic reactions.

DISPOSAL

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

Storage:

Processing:

at 23 °C ± 2 °C

Manufacturing process

Data processing and creation of the support structure according to the specifications of the CAD software manufacturer

Building process

Creation of a print job in compliance with the machine and material parameters.

Finishing process

After raising the platform, a dripping time of approx. 10 minutes is recommended. If possible, the finishing should take place immediately after the building process.

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*The set of parameters includes all relevant material- and printer specific information

CLEANING EQUIPMENT

Cleaning unit Manufacturer, Model	Cleaning process
Rapidshape RS Wash	Use the following settings: DETAX Freeprint-tryin 385 Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.
Ultrasonic bath Bandelin Sonorex	Clean the parts with isopropyl alcohol (purity \geq 98 %) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air. The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity \geq 98 %) for 3 minutes. Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.

LIGHT CURING EQUIPMENT

Light Curing unit Manufacturer, Model	Curing process
Dentalfarm Photopool	2 x 3 min, progressive + N ₂ , turnaround components after 3 min
NK Optik Otoflash G171	2 x 2000 flashes under inert gas, turnaround components after 2000 flashes
NK Optik Otoflash 250/500	4000 flashes under inert gas @ 15 Hz
Rapidshape RS cure	Use the following settings: DETAX Freeprint-tryin 385

DETX

DETAX GmbH
Carl-Zeiss-Str. 4 · 76275 Ettlingen/Germany
☎ +49 7243/510-0 · 📧 +49 7243/510-100
www.detax.com · post@detax.com

Made in
Germany
09/2023
3.1