blue eco scan

DE

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EN

Instructions for use

Scannable laboratory putty, addition curing, smooth, non-sticky consistency

For the powderless 3D data acquisition of the antagonists for the CAD/CAM/CIM technologies (e.g. strip light projection, laser triangulation and CEREC procedures), excellent resolution for artifactfree image reproduction with high brightness and strong contrast. Easy mixing due to contrast colours and smooth consistency. Can be cut and milled very well, no breaking out of thin contact point areas.

Mixing and dispensing

Remove components A (white) and B (light grey) with the corresponding measuring scoop and mix quickly in the ratio 1:1 (approx. 30 seconds), until a homogeneous mixture or uniform colour is achieved.


Homogeneous kneading

Cutting

Preparing for scan

Screenshot

Important working hints

- For further work on the registration, observe the instructions of the scan system manufacturer.
- Clean and dry the registration before scanning.
- Only for dental **extraoral** use by qualified personnel.
- Make sure that only base and catalyst components from the same batch number are used.
- Do not bring the components A + B into contact before mixing, do not interchange measuring scoops and jar lids.
- Do not remove material by hand, do not wear latex gloves.
- Please follow the instructions of the safety data sheet!

1 min. 30 sec.

Working time

4 min. 30 sec.

Setting time

Indications

for use:

- Occlusal registration for powderless 3D data acquisition
- Fabrication of small scannable model segments
- Conventional indications for laboratory putty

Technical Data:

- **Mixed volume:** 2 x 800 g (jars)
- **Mixing ratio:** 1:1
- **Colour code:** base: light grey catalyst: white
- **Mixing time:** approx. 30 sec.
- **Working time:** approx. 1 min. 30 sec.*
- **Setting time:** approx. 4 min. 30 sec.*
- **Final hardness:** approx. 85 Shore A
- **Linear dimensional change:** approx. 0.2 %
- **Application:** At 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity
- **Stockage:**

* from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard a. m. times.

Ordering

information:

blue eco scan 03118

Standard packing
base, jar of 800 g
catalyst, jar of 800 g
2 measuring scoops

blue eco 02038

Standard packing
base, jar of 800 g
catalyst, jar of 800 g
2 measuring scoops

Mega pack 02467

base, bucket of 5 kg
catalyst, bucket of 5 kg
2 measuring scoops

blue eco stone 03117

Standard packing
base, jar of 800 g
catalyst, jar of 800 g
2 measuring scoops

Eco pack 02517

base, jar of 1400 g
catalyst, jar of 1400 g
2 measuring scoops

Mega pack 02985

base, bucket of 5 kg
catalyst, bucket of 5 kg
2 measuring scoops

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