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EN

Instructions for use

Light-curing, mouldable resin that burns out without residue

Light-curing, ready-to-use, single component material for use as a modelling gel and paste for multiple indications in dental engineering, especially mould casting. The gel is applied straight from the metering syringe and the paste can easily be moulded into shape thanks to its kneadable consistency. **EASYFORM LC** can be polymerised in all currently available light-curing devices. After the final curing the material can easily be trimmed both by the micro engine and by milling devices. **EASYFORM LC** burns out without residue during the preheating of the casting mould.

Fields of application

Fixed dentures

- Modelations for casting of crowns and bridges, inlays and onlays.

Combination technique

- Primary and secondary components for the telescopic and tapered crown technique, customised attachments, secondary components for groove-shoulder-pin attachments, stress distribution arms.

Implantology

- Models for casting implant abutments, implant bars and implant superstructures and tertiary frameworks. Fabrication of transfer and fixation keys and splinting or fixation of customised parts of implants or segments

Joining technique and fixations

- Fixation of parts or segments by soldering or welding
- Fixation of prefabricated attachments, general splinting or fixation for repairs or trying in dentures

Transfer and bite check

- Transfers splints and keys, that is insertion guides, occlusal bite splints, stops and check bites

1. Preparation

Insulate the surfaces of the plaster with standard alginates. Vaseline is suitable for insulating polyurethane modelling materials. As a rule, no further insulation is required on polished metal or ceramic surfaces. If necessary, apply a wafer thin layer of Vaseline or oil, blowing off any excess oil.

2. Application

Remove the protective cap from the application syringe. Fit black application needles onto the application syringe and lock in place. Press the piston lightly to dispense and model the material. The material flows easily and can be guided while the tip is being moved in the material but stops flowing as soon as the tip is stationary. When modelling parts with thicker walls, apply the material in several thinner layers of approx. 1 mm. Before applying the next layer an intermediate polymerisation of 10 sec. duration is recommended for the handheld lamp or 60 sec. for the laboratory apparatus. This prevents incomplete curing, excessive release of heat or stresses. After use reseal the syringe with the protective cap. Remove the paste with a spatula or dental instrument. Always replace the lid of the paste tin immediately after removing the material in order to protect the contents from the light. Preparatory shaping can be carried out by hand. Use the paste for applying larger amounts of the material quickly and effectively and the gel for corrections and jointing. After the final polymerisation the model can be trimmed or shaped in the milling device. The model should be covered with a thin layer of modelling wax for large scale work.

3. Curing

EASYFORM LC gelpaste can be polymerised in all currently available light-curing devices with a light spectrum of 320 nm to 435 nm. Devices with a curing chamber equipped with a cooling system are also suitable. The average polymerisation times in light-curing devices with UV/UV-A lamps are approx. 3 min. 30 sec. to 5 min., for halogen lamps approx. 2 min., for stroboscopic lamps approx. 1 min. 30 sec. and for LED lamps approx. 20 sec. It should be noted that the polymerisation time is always dependent on the lamp in the device. Consistent curing times can be guaranteed only by the use of regularly maintained devices. Follow the manufacturer's instructions.

Important points to observe when using the material

- In order to achieve optimum casting results and cast surfaces, the oxygen inhibition layer must be completely removed before model is immersed in isopropanol.
- Modelling of **EASYFORM LC paste** is greatly facilitated by coatings fingers or instruments with petroleum jelly.
- The occasional forming of separations on the paste surface in the tin is due to changes in temperature and does not represent an impairment of the paste.
- Trimming and polishing should be carried out only at workplaces equipped with dust extraction units.
- Avoid contact of the product with the mucous membranes and eyes. Wear suitable eye protection and gloves! If accidental contact occurs, rinse under running water and seek medical assistance if necessary. Material that has not been cured can cause irritation and sensitisation if it comes into contact with the skin and mucous membranes.
- Carefully reseal tin and syringe immediately after use.
- The product has been developed for use in dental applications and must be used only for its intended purpose by qualified personnel.
- Please follow the instructions of the safety data sheet!

Technical data:

- **Volumes**
3 g syringe/30 g tin
- **Dispensing and application**
Direct application by application needles or spatula
- **Product colour**
transparent red
- **Curing**
Light-curing devices with a light spectrum of 320 nm to 435 nm
- **UV/UV-A lamps**
approx. 3 min. 30 sec. to 5 min.
- **Halogen lamps**
approx. 2 min.
- **Stroboscopic lamps**
approx. 1 min. 30 sec.
- **LED lamps**
approx. 20 sec.
- **Intermediate polymerisation per 1 mm layer**
Hand held lamp approx. 10 sec.
Laboratory apparatus approx. 60 sec.
- **Application**
23 °C ± 2 °C, 50 ± 5% relative humidity
- **Storage**

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